



## Specifications:

- Patient weight capacity up to 1000 pounds
- Rugged reinforced frame and litter top
- Three position vertical side rails
- Foot pedal control hydraulic fowler back rest
- Four wheel central locking
- Full perimeter non-marring bumper
- Four IV receptacles
- 4" hinged high density mattress pad
- Litter top usable patient surface measures 37" x 80"
- Overall bumper to bumper width 40" for doorway access
- Litter top to floor height measures 29-3/8"
- Beachwood powder coat finish
- 8" All surface casters
- Shipping weight: 490 Lbs.

## Options:

- Utility Shelf
- IV Rod Holder
- Oxygen Holder
- Restraint Straps
- Padded Arm Board
- Chart Holder