



Y O U T H B E D



Specifications

Manual: 3 Manual Cranks: Head Gatch, Knee Gatch, and Hi-Lo Positioning

Electric: Hand Control for the above positions (UL Approved)

4 IV Hook Holes

5" Braking Casters

Mattress Stop

Solid Pan Sleeping Surface (non-helical design)

Optional 5" Hinged Comfort Foam Mattress (Soft, Blue Cover)

Dimensions

Overall Length: 72.5"

Overall Width: 38" with Rails UP

Overall Height at HI position: 41.5"

Overall Height at LO position: 32"

Rail Dimension (Height): 25.75"

Sleeping Surface 33.5" x 64.5"

Mattress Surface Ht at HI position: 26.75"

Mattress Surface Ht at LO position: 17.25 (no Mattress)

For more information, either call **(800) 274-2742** or email us at **info@nkmedicalproducts.com**



NK Medical Products, Inc.; 10123 Main St, PO Box 627, Clarence, NY 14031

<http://www.nkmedicalproducts.com>